



Recombinant Human Putative SCAN domain-containing protein 2 (SCAND2)

Product Code	CSB-EP861977HU
Abbreviation	SCAND2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZW5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAVAVDQQIQ TPSVQDLQIV KLEEDSHWEQ EISLQGNYPG PETSQCQSFWH FRYQEASRPR EALLQLQKLC CQWLRPEKCT KEQILELLVL EQFPTVLLQE IQIWVRQQHP ESGEEAVALV EDLQKEPGRQ RLEPRARPSG RTPPAQLRSP WPMTAAGPAS RARASETGST ASCAGRWRTC CAAAAAPSAA RSASARTGRS TSSCARAARA PSATEGALTR TPAPRRPLQR RRPGTGPWRP GRQRGAGTAP PGTQPRQRPS PRPTPRRPRP RLARPRAGQK PWLLVIRSME RVMYVMLIIQ METEDV
Source	E.coli
Target Names	SCAND2P
Protein Names	Recommended name: Putative SCAN domain-containing protein 2
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.