



Recombinant Human UPF0696 protein C11orf68 (C11orf68)

Product Code	CSB-EP862008HU-B
Abbreviation	C11orf68
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3H3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEPGEELEEE GSPGGREDGF TAEHLAAEAM AADMDPWLVF DARTTPATEL DAWLAKYPPS QVTRYGDPGS PNSEPVGWIA VYGQGYSPNS GDVQGLQAAW EALQTSGRPI TPGTLRQLAI THHVLSGKWL MHLAPGFKLD HAWAGIARAV VEGQLQVAKV SPRAKEGGRQ VICVYTDDFT DRLGVLEADS AIRAAGIKCL LTYKPDVYTY LGIYRANRWH LCPTLYESRF QLGGSSARGSR VLDRANNVEL T
Source	E.coli
Target Names	C11orf68
Protein Names	Recommended name: UPF0696 protein C11orf68 Alternative name(s): Basophilic leukemia-expressed protein Bles03 Protein p5326
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.