



Recombinant Human Tyrosine-protein kinase Srms (SRMS)

Product Code	CSB-EP862011HU-B
Abbreviation	SRMS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3Y6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEPFLRRRLA FLSFFWDKIW PAGGEPDHGT PGSLDPNTDP VPTLPAEPCS PFPQLFLALY DFTARCGGEL SVRRGDRLCA LEEGGYIFA RRLSGQPSAG LVPITHVAKA SPETLSDQPW YFSGVSRTQA QQLLLSPPNE PGAFLIRPSE SSLGGYSLSV RAQAKVCHYR VSMAADGSLY LQKGR LFPGL EELLTYKAN WKLIQNPLLQ PCMPQKAPRQ DVWERPHSEF ALGRKLGEGY FGEVWEGLWL GSLPVAIKVI KSANMKLTDL AKEIQTLKGL RHERLIRLHA VCSGGEPVYI VTELMRKGNL QAFLGTPEGR ALRLPPLLGF ACQVAEGMSY LEEQRVVHRD LAARNVLVDD GLACKVADFG LARLLKDDIY SPSSSSKIPV KWTAPEAANY RVFSQKSDVW SFGVLLHEVF TYGQCPYEGM TNHETLQQIM RGYRLPRPAA CPAEVYVLMML ECWRSSPEER PSFATLREKL HAIHRCHP
Source	E.coli
Target Names	SRMS
Protein Names	Recommended name: Tyrosine-protein kinase Srms EC= 2.7.10.2
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.