



# Recombinant Human Zinc finger and SCAN domain-containing protein 16 (ZSCAN16)

<b>Product Code</b>	CSB-EP862018HU-B
<b>Abbreviation</b>	ZSCAN16
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H4T2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTALEPEDQ KGLLIKAED HYWGQDSSSQ KCSPHRRELY RQHFRKLCYQ DAPGPREALT QLWELCRQWL RPECHTKEQI LDLLVLEQFL SILPKDLQAW VRAHHPETGE EAVTVLEDLE RELDEPGKQV PGNSERRDIL MDKLAPLGRP YESLTVQLHP KKTQLEQEAG KPQRNGDKTR TKNEELFQKE DMPKDKEFLG EINDRLNKDT PQHPKSKDII ENGRSEWQQ RERRRYKCDE CGKSFSHSSD LSKHRRTHTG EKPYKCDECG KAFIQRSHLI GHHRVHTGVK PYKCKECKGD FSGRTGLIQH QRIHTGEKPY ECDECGRPFR VSSALIRHQR IHTANKLY
<b>Source</b>	E.coli
<b>Target Names</b>	ZSCAN16
<b>Protein Names</b>	Recommended name: Zinc finger and SCAN domain-containing protein 16 Alternative name(s): Zinc finger protein 392 Zinc finger protein 435
<b>Expression Region</b>	1-348
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.