



# Recombinant Human Zinc finger and SCAN domain-containing protein 16 (ZSCAN16)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP862018HU   |
| <b>Abbreviation</b>      | ZSCAN16  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q9H4T2   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MTTALEPEDQ KGLLIKAED HYWGQDSSSQ KCSPHRRELY RQHFRKLCYQ<br>DAPGPREALT QLWELCRQWL RPECHTKEQI LDLLVLEQFL SILPKDLQAW<br>VRAHHPETGE EAVTVLEDLE RELDEPGKQV PGNSERRDIL MDKLAPLGRP<br>YESLTVQLHP KKTQLEQEAG KPQRNGDKTR TKNEELFQKE DMPKDKEFLG<br>EINDRLNKDT PQHPKSKDII ENGRSEWQQ RERRRYKCDE CGKSFSHSSD<br>LSKHRRTHTG EKPYKCDECG KAFIQRSHLI GHHRVHTGVK PYKCKECKGD<br>FSGRTGLIQH QRIHTGEKPY ECDECGRPFR VSSALIRHQR IHTANKLY |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | ZSCAN16  |
| <b>Protein Names</b>     | Recommended name: Zinc finger and SCAN domain-containing protein 16<br>Alternative name(s): Zinc finger protein 392 Zinc finger protein 435  |
| <b>Expression Region</b> | 1-348  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.    |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |