



Recombinant Human Coiled-coil domain-containing protein 86 (CCDC86)

Product Code	CSB-EP862027HU-B
Abbreviation	CCDC86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6F5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDTPLRRSRR LGGLRPESPE SLTSVSRTRR ALVEFESNPE ETREPGSPPS VQRAGLGSPE RPPKTSPGSP RLQQGAGLES PQQQPEPGAA SPQRQQDLHL ESPQRQPEYS PESPRCQPKP SEEAPKCSQD QGVLASELAQ NKEELTPGAP QHQLPPVPGS PEPYPGQQAP GPEPSQPLLE LTPRAPGSPR GQHEPSKPPP AGETVTGGFG AKKRKGSSSQ APASKKLNKE ELPVIPKGKP KSGRVWKDRS KKRFSQMLQD KPLRTSWQRK MKERQERKLA KDFARHLEEE KERRRQEKKQ RRAENLKRRL ENERKAEVVQ VIRNPAKLKR AKKKQLRSIE KRDTLALLQK QPPQQPAAKI
Source	E.coli
Target Names	CCDC86
Protein Names	Recommended name: Coiled-coil domain-containing protein 86 Alternative name(s): Cytokine-induced protein with coiled-coil domain
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.