



# Recombinant Human RNA polymerase II-associated protein 3 (RPAP3)

<b>Product Code</b>	CSB-BP862030HU
<b>Abbreviation</b>	RPAP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H6T3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TSANKAIEL QLQVKQNAEE LQDFMRDLEN WEKDIKQKDM ELRRQNGVPE ENLPPIRNGN FRKKKKGKAK ESSKKTREEN TKNRIKSYDY EAWAKLDVDR ILDELKDDDS THESLSQESE SEEDGIHVDS QKALVLKEKG NKYFKQGKYD EAIDCYTKGM DADPYNPVLV TNRASAYFRL KKFVAESDC NLAVALNRSY TKAYSRRGAA RFALQKLEEA KKDYERVLEL EPNNFEATNE LRKISQALAS KENSYPKEAD IVIKSTEGER KQIEAQQNKQ QAISEKDRGN GFFKEGKYER AIECYTRGIA ADGANALLPA NRAMAYLKIY KYEEAEKDCT QAILLDGSYS KAFARRGTAR TFLGKLNEAK QDFETVLLLE PGNKQAVTEL SKIKKELIEK GHWDDVFLDS TQRQNVVKPI DNPPHPGSTK PLKKVIIIEET GNLIQTIDVP DSTTAAAPEN NPINLANVIA ATGTTSSKNS SQDDLFPSTSD TPRAKVLKIE EVSDTSSLQP QASLKQDVCQ SYSEKMPIEI EQKPAQFATT VLPPIPANSF QLESDFRQLK SSPDMLYQYL KQIEPSLYPK LFQKNLDPDV FNQIVKILHD FYIEKEKPLL IFEILQRLSE LKRFDMAVMF MSETEKKIAR ALFNHIDKSG LKDSSVEELK KRYGG
<b>Source</b>	Baculovirus
<b>Target Names</b>	RPAP3
<b>Protein Names</b>	Recommended name: RNA polymerase II-associated protein 3
<b>Expression Region</b>	2-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes an RNA polymerase II-associated protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.