



Recombinant Human Methyltransferase-like protein 7A (METTL7A)

Product Code	CSB-BP862041HU
Abbreviation	METTL7A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8H3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	W ICKKWFPYFL VRFTVIYNEQ MASKKRELFS NLQEFAGPSG KLSLLEVGCG TGANFKFYPP GCRVTCIDPN PNFEKFLIKS IAENRHLQFE RFVVAAGENM HQVADGSVDV VVCTLVLCVS KNQERILREV CRVLRPGGAF YFMEHVAEAC STWNYFWQQV LDPAWHLLFD GCNLTRESWK ALERASFSKL KLQHIQAPLS WELVRPHIYG YAVK
Source	Baculovirus
Target Names	METTL7A
Protein Names	Recommended name: Methyltransferase-like protein 7A EC= 2.1.1.- Alternative name(s): Protein AAM-B
Expression Region	30-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.