



Recombinant Human Plexin-A4 (PLXNA4), partial

Product Code	CSB-EP862071HU-B
Abbreviation	PLXNA4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCM2
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LLTRQPAPLSQKQRSFVTRGEPAEAGFNHLVVDERTGHIYLGAVNRIYKLSSDL KVLVTHETGPDENPKCYPPRIVQTCNEPLTTTNNVNMKMLLIDYKENRLIACGS LYQGICKLLRLEDLFLKLGEPYHKKEHYLSGVNESGVSFVGVIVSYSNLDDKLFAT AVDGKPEYFPTISSRKLTKNSEADGMFAYVFHDEFVASMIPSDTFTIIPDFDIY YVYGFSSGNFVYFLTLQPEMVSPPGSTTKEQVYTSKLVRLCKEDTAFNSYVEV PIGCERSGVEYRLLQAAYLSKAGAVLGRTLGVHPDDDLLFTVFSKGQKRKMKS LDESALCIFILKQINDRIKERLQSCYRGEGLDLAWLKVKDIPCSSALLTIDNFC GLDMNAPLGVSDMVRGIPVFTEDRDRMTSVIAYVYKNHSLAFVGTKSGKLLKKS FGTGPQGGITQEWIGVEGDPGGANIASQEQLCVYLQCSSHKAISDQRVQPLL CCFLNVPGNS
Source	E.coli
Target Names	PLXNA4
Protein Names	Recommended name: Plexin-A4
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial of Isoform 3
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.