



Recombinant *Leifsonia aquatica* Levodione reductase (lvr)

Product Code	CSB-EP862753LEAR
Abbreviation	lvr
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LBG2
Product Type	Recombinant Protein
Immunogen Species	<i>Leifsonia aquatica</i> (<i>Corynebacterium aquaticum</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTATSSPTTR FTDRVVLITG GGSGLGRATA VRLAAEGAKL SLVDVSSEGL EASKAAVLET APDAEVLTTV ADVSDEAQVE AYVTATTERF GRIDGFFNNA GIEGKQNPT E SFTAAEFDKV VSINLRGVFL GLEKVLKIMR EQGSGMVVNT ASVGGIRGIG NQSGYAAAKH GVVGLTRNSA VEYGRYGIRI NAIAPGAIWT PMVENS MKQL DPENPRKAAE EFIQVNPSKR YGEAPEIAAV VAFLLSDDAS YVNATVVPID GGQSAAY
Source	E.coli
Target Names	lvr
Protein Names	Recommended name: Levodione reductase EC= 1.1.1.- Alternative name(s): (6R)-2,2,6-trimethyl-1,4-cyclohexanedione reductase
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.