



Recombinant Beta-1,3-xylanase TXYA (txyA)

Product Code	CSB-MP862754VFF
Abbreviation	txyA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LCB9
Product Type	Recombinant Protein
Immunogen Species	Vibrio sp.
Purity	>85% (SDS-PAGE)
Sequence	LDGKLVPD QGILVSVGQD VDSVNDYSSA MGTTPAGVTN YVGIVNLDGL STDADAGAGR NNIVELANQY PTSALIVGVS MNGEVQNVAN GQYNANIDTL IRTLGEFDRP VYLRWAYEVD GPWNGHNTED LKQSFHRVYQ RIRELGYADN ISMVWQVASY CPTAPGQLGT WWPGDVVVDW VGLSYFAPQD CNWDRVNEAA QWARSHNKPL FINESPQRY QLADLTYSTD PAKGTNRQAK TDQQIWSEWF EPPFQFMVDN QDILKGFTYI NADWDSQWRW AAPYNEGYWG DSRVQVIPYI KQKWQETLSD PKFIRHSEDL FAQLGYGNSD GGNGGDNGGD NGGDNGGETP ENCTDDFNFN YVSDNEIEVY HVDKGWSAGW NYLCLDDYCL SGTKSNGAFS RSFSAQLGQT YKMTFKVEDI TGQQQQIIDK TVTFTNQVCN
Source	Mammalian cell
Target Names	txyA
Protein Names	Recommended name: Beta-1,3-xylanase TXYA EC= 3.2.1.32 Alternative name(s): Endo-1,3-beta-xylanase
Expression Region	23-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.