



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g55890, mitochondrial (At1g55890)

Product Code	CSB-MP862782DOA
Abbreviation	At1g55890
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LG23
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	F SSAATVVSEP TAVTAAISPP QKSLTSLVNG ERNPKRIVEK FKKACESERF RTNIAVYDRT VRRLVAAKRL HYVEEILEEQ KKYRDMSEK FAARIISLYG KAGMFENAQK VFEEMPNRDC KRSVLSFNAL LSAYRLSKKF DVVEELFNEL PGKLSIKPDI VSYNTLIKAL CEKDSLPEAV ALLDEIENKG LKPDIVTFNT LLLSSYLKGG FELGEEIWAK MVEKNVAIDI RTYNARLLGL ANEAKSKELV NLFGELKASG LKPDVFSFNA MIRGSINEGK MDEAEAWYKE IVKHGYRDPK ATFALLLPAM CKAGDFESAI ELFKETFSKR YLVGQTTLLQQ LVDELVKGSK REEAEEIVKI AKTNDFLKLL LNLPSQEE
Source	Mammalian cell
Target Names	At1g55890
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At1g55890, mitochondrial
Expression Region	20-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.