



# Recombinant Arabidopsis thaliana

## Uncharacterized protein At3g28850 (At3g28850)

<b>Product Code</b>	CSB-EP862786DOA-B
<b>Abbreviation</b>	At3g28850
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LH89
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGCASSKHRK RCLHCRRGYS PPVDVQRSHS VHHASQNSED SCHMVALSSS SLGSLKLCDS SFGHNHKKHLA DFSEKLVSGE SVKTGNGFGP NVVREKSDKE KSNLELQAKL MEAKVWSSMM NEKIPKIVPK TPVTPPGEP ETINTWEMMD GLEDVLSPLR SPNHVKSFSF DVGPNNGGKSN GSVKPVWLQM EEEEEGFEDF DPEIISFRK SLQELPSDHP FHISNHDFEL KPRFNFSDEE KEEEEEQSVGK ERVILYFTSL RGIRKTYEES CDVRVILKSL GIRVDERDVS MHSGFKDELK ELLGEKFNKG VGITLPRVFL GRKYIGGAEE IRKLNEDGKL EKLLGGCERV EENQNGNGLE CEACGDVRFV PCETCSGSCK VYYEYEDDDD DDDEGDDDES VKEEREYGFQ TCPDCNENGL IRCPVCCD
<b>Source</b>	E.coli
<b>Target Names</b>	At3g28850
<b>Protein Names</b>	Recommended name: Uncharacterized protein At3g28850
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.