



Recombinant Arabidopsis thaliana

Uncharacterized protein At3g28850 (At3g28850)

Product Code	CSB-MP862786DOA
Abbreviation	At3g28850
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LH89
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGCASSKHRK RCLHCRRGYS PPVDVQRSHS VHHASQNSD SCHMVALSSS SLGSLKLCDS SFGHNHKKHLA DFSEKLVSGE SVKTGNGFGP NVVREKSDKE KSNLELQAKL MEAKVWSSMM NEKIPKIVPK TPVTPPGEP ETINTWEMMD GLEDVLSPLR SPNHVKSFSF DVGPNNGGKSN GSVKPVWLQM EEEEEGFEDF DPEIISFRK SLQELPSDHP FHISNHDFEL KPRFNFSDEE KEEEEEQSVGK ERVILYFTSL RGIRKTYEES CDVRVILKSL GIRVDERDVS MHSGFKDELK ELLGEKFNKG VGITLPRVFL GRKYIGGAEE IRKLNEDGKL EKLLGGCERV EENQNGNGLE CEACGDVRFV PCETCSGSCK VYYEYEDDDD DDDEGDDDES VKEEREYGFQ TCPDCNENGL IRCPVCCD
Source	Mammalian cell
Target Names	At3g28850
Protein Names	Recommended name: Uncharacterized protein At3g28850
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.