



Recombinant Arabidopsis thaliana Peroxidase 32 (PER32)

Product Code	CSB-YP862787DOA
Abbreviation	PER32
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LHB9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	Q LTPTFYDNTC PSVFTIVRDT IVNELRSDPR IAASILRLHF HDCFVNGCDA SILLDNTTSF RTEKDAAPNA NSARGFPVID RMKAAVETAC PRTVSCADIL TIAAQQAVNL AGGPSWRVPL GRRDSLQAFF ALANTNLPAP FFTLPQLKAS FQNVGLDRPS DLVALSGGHT FGKNQCQFIM DRLYNFSNTG LPDPTLNNTY LQTLRGQCPR NGNQTVLVDF DLRTPTVFDN KYYVNLKELK GLIQTQELF SSPNATDTIP LVREYADGTQ KFFNAFVEAM NRMGNITPLT GTQGQIRQNC RVVNSNSLLH DVVEIVDFVS SM
Source	Yeast
Target Names	PER32
Protein Names	Recommended name: Peroxidase 32 Short name= Atperox P32 EC= 1.11.1.7 Alternative name(s): ATP16a PRXR3
Expression Region	30-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.