



Recombinant Arabidopsis thaliana Putative F-box protein At3g23970 (At3g23970)

Product Code	CSB-EP862797DOA
Abbreviation	At3g23970
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LIR0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNIPPELTFE VLVRLPLKSL ARFRSVRKEW KLVIDSEFFR DCFMSHNSSS VSWSIIQTRP HKLTLEIVGH HGCKTWGLTR SPGSLVGFFA DTTIKRLRVL ACTDGLVLIY TETSYGTPMH YVGNPLFQEW CPIPLPSYFY LQSVERIRNH QRFNDSALVI KMQSGAVVSY KVVWLFPRNS TTIDFLIYSS DTGMWEARIA TCLHSSFWFS HHKPIALNGI LHWLCNFSTS FVAYDFYGGH DDYACHIISF PDCEKDDDQL RFRRTLSTS DGSIVYFNEF GENGNRRLRV WRLVKYTGGP EAWQLFWEVS LASVTKLGID YFPVVMHPLK SEIYLWSRN KKGMLVFNLR TNVFLHKES EDETKCMDGC TFSFNRCSEY METIHNYGGP NYLLASQYVL PRWLHRLPRP QPS
Source	E.coli
Target Names	At3g23970
Protein Names	Recommended name: Putative F-box protein At3g23970
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.