



Recombinant Arabidopsis thaliana Putative SNAP25 homologous protein SNAP30 (SNAP30)

Product Code	CSB-MP862819DOA
Abbreviation	SNAP30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LMG8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MFGFFKSPGN NKLPNESSNN KGGTITAGRR TSSEPILITP DFDDDDKYKN GFNDSSGGLQS QTTEELEKYA VYKAEETTKG VNNCLKIAED IRSDGARTLE MLHQQGEQIN RTHEMAVDMD KDLSRGEKLL NNLGGMFSKP WPKPKTKNIT GPMITPDKPS KKSENHKEER EKLGLGAKGR SSSQPALDQP TNALQKVEQE KAKQDDGLSD LSDILGDLKS MAVDMGSEID KQNKALDHLG DDVDELNSRV QGANQRARHL LSK
Source	Mammalian cell
Target Names	SNAP30
Protein Names	Recommended name: Putative SNAP25 homologous protein SNAP30 Short name= AtSNAP30 Alternative name(s): Synaptosomal-associated protein SNAP25-like 3
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.