



Recombinant Arabidopsis thaliana Transcription factor TCP23 (TCP23)

Product Code	CSB-MP862845DOA
Abbreviation	TCP23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LQF0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MESHNNNQSN NNTTGS AHLV PSMGPISGSV SLTTTAPNST TTTVTA AKTP AKRPSKDRHI KVDGRGRRIR MPAICAAARVF QLTRELQHKS DGETIEWLLQ QAEPAAIAAT GTGTIPANIS TLNISRSSG STLSAPLSKS FHMGRAAQNA AVFGFQQQLY HPHHITDSS SSSLPKTFRE EDLFKDPNFL DQEPGSRSPK PGSEAPDQDP GSTRSRTQNM IPPMWALAPT PASTNGGSAF WMLPVGGGGG PANVQDPSQH MWAFNPGHYP GRIGSVQLGS MLVGGQQLGL GVAENNNLGL FSGGGGDGGR VGLGMSLEQK PQHQVSDHAT RDQNPTIDGS P
Source	Mammalian cell
Target Names	TCP23
Protein Names	Recommended name: Transcription factor TCP23
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.