



Recombinant Arabidopsis thaliana Peroxisomal (S)-2-hydroxy-acid oxidase GLO1 (GLO1)

Product Code	CSB-MP862860DOA
Abbreviation	GLO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LRR9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEITNVTEYD AIAKQKLPKM VYDYYASGAE DQWTLQENRN AFARILFRPR ILIDVSKIDM TTTVLGFKIS MPIMVAPTAM QKMAHPDGEY ATARAASAAG TIMTLSSWAT SSVEEVASTG PGIRFFQLYV YKNRNVVEQL VVRAERAGFK AIALTVDTPR LGRRES DIKN RFTLPPNLT KNFEGLDLGK MDEANDSGLA SYVAGQIDRT LSWKDVQWLQ TITKLPIVK GVLGTGEDARI AIQAGAAGII VSNHGARQLD YVPATISALE EVVKATQGRI PVFLDGGVRR GTDVF KALAL GASGIFIGRP VVFSLAAEGE AGVRKVLQML RDEFELTMAL SGCRSLKEIS RNHITTEWDT PRPSARL
Source	Mammalian cell
Target Names	GLO1
Protein Names	Recommended name: Peroxisomal (S)-2-hydroxy-acid oxidase GLO1 EC= 1.1.3.15 Alternative name(s): Glycolate oxidase 1 Short name= AtGLO1 Short name= GOX 1 Short chain alpha-hydroxy acid oxidase GLO1
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.