



Recombinant Arabidopsis thaliana Putative F-box protein At3g25750 (At3g25750)

Product Code	CSB-BP862863DOA
Abbreviation	At3g25750
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LS04
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEKTEWSDLP EELLDLIANR YSSNIDVLRI RSTCKSWRSA VAMSKERLQF RFERYLPTSN KKIKAHLSPT TFFRITLPSS CPNKGWLVRT RQASKMYRKI TLLCPLSGER ITRSHQTLDL LKVGVSSEIRQ SYEIQIFDGL KDEKIPLDSE IFSNIYKNSD KIPSVNYHDN KIWCCKTREG SRSWTKIKNQ VEDFSDIILH MGRIYAVDLK GAIWWISLSQ LTIVQQTSST PLDYYKYDSC QDTRLVEYCG DLCIVHELSE TRNHIQRTVG FKVKYKMDDEL AKWVEVSCLG DNTLIVACNS CFTVVASEYH GCLKNSIYFS YYDVKKAENI KVFKLDDGSI TQRTDISSQS CFHMFSLPFL NY
Source	Baculovirus
Target Names	At3g25750
Protein Names	Recommended name: Putative F-box protein At3g25750
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.