



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At3g17540 (At3g17540)

Product Code	CSB-BP862896DOA
Abbreviation	At3g17540
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUP4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRKTMVISDL PHEIESEILS RVPTKSLAKL HTTCKRWYAL FRDPRFVKKN FGKSERRML HSNFGVYKIT DDLHGILNSG DPSLEFTSKL SNLKISEDLT ITKIFHCDGL ILCSTKENTR LVVWNPCTGQ TRWIKPSKRY RSDDSYCLGY VNSKSSYHNY KILRYCFYYN DQDACVSEFE IYDFSSSWR VLDDYCTREW GLFCHGMSLK GNTYFVAGEK ETGFFMLYFD FKTERFERLP LPYQSFDS ED TAVLSIVGGE KLAVLHQNIQ SFSNEMRIWV TNKIDEAKDL TWSNFFLA VD YDIFNLPSVN NVTSFLLDEE NKVAVCCDRH IDDEDKTRIIY IVGVDLYKEV YKERTKGAHF NWPLLSYLP SLVHIQEKYP KDKRKGKRR LVKRIT
Source	Baculovirus
Target Names	At3g17540
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At3g17540
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.