



Recombinant Arabidopsis thaliana 60S acidic ribosomal protein P3-2 (RPP3B)

Product Code	CSB-MP862904DOA
Abbreviation	RPP3B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LVC9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGVFSFVCKS KGGWETAKQH EGDLEASASS TYDLQRKLVQ TALSADSSGG VQSSFSLVSP TSAVFVVVIG GGGGGGFAAG GGAAAGGGGG GEAAAATKEE EKKKEESEEE EGDFGFDLFG
Source	Mammalian cell
Target Names	RPP3B
Protein Names	Recommended name: 60S acidic ribosomal protein P3-2
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.