



# Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At3g52680 (At3g52680)

<b>Product Code</b>	CSB-BP862922DOA
<b>Abbreviation</b>	At3g52680
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LXJ6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MISPTMEQGE SLRIRVVMGK DRISLPDGL LKILSSLPT NIVVATSVLS KQWRSLWKLV PNLEFDSDDY ESEHYTFSEI VCKSFLSHKA PVLESFRLKF VNFNPVDIGL WVGIAFSRHL RELVLDYFPA ELGKGVTFTF PSSLCTCNTL ETLKLVLCIL VDIPSPVLMK SLRTLHLEFV RYKDESSVRN LLSGCPGLEE LRLYRGDDSD IKVFTIEVPS LQRLTIHDNN DGPEFWGYVI NAPFLKYLLI EELRCPEFCL NAPELVEANI AEVTSITIEK FLGSFTSVSR LLLNLSPLKI TYPTGSMFYQ LVSLEMYTRE AEWWNLLTLM LENS PKLQVL KLTDRSQNFH KDGLVSGKWN EPKDVPECLL SQLETFVWRR FDWGREEEKE IATYILKNGR RLKKATFSTN PIESEELNKL KERRKVLNGL DGVVRASNCSQLVFKFDPSY LGSDSP
<b>Source</b>	Baculovirus
<b>Target Names</b>	At3g52680
<b>Protein Names</b>	Recommended name: F-box/FBD/LRR-repeat protein At3g52680
<b>Expression Region</b>	1-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.