



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At3g52680 (At3g52680)

Product Code	CSB-MP862922DOA
Abbreviation	At3g52680
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXJ6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MISPTMEQGE SLRIRVVMGK DRISLPDGL LKILSSLPT NIVVATSVLS KQWRSLWKL V PNFVDSDDY ESEHYTFSEI VCKSFLSHKA PVLESFRLKF VNFNPVDIGL WVGIAFSRHL RELVLDYFPA ELGKGVTF TF PSSLCTCNTL ETLKLVL CIL VDIPSPVLMK SLRTLHLEFV RYKDESSVRN LLSGCPGLEE LRLYRGDDSD IKVFTIEVPS LQRLTIHDNN DGPEFWGYVI NAPFLKYLLI EELRCPEFCL NAPELVEANI AEVTSITIEK FLGSFTSVSR LLLNLSPLKI TYPTGSMFYQ LVSLEMYTRE AEWWNLLTLM LENS PKLQVL KLTDRSQNFH KDGLVSGKWN EPKDVPECLL SQLETFVWRR FDWGREEEKE IATYILKNGR RLKKATFSTN PIESEELNKL KERRKVLNGL DGVVRASN SC QLVFKFDPSY LGSDSP
Source	Mammalian cell
Target Names	At3g52680
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At3g52680
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.