



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At3g44090 (At3g44090)

Product Code	CSB-MP862923DOA
Abbreviation	At3g44090
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXQ3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNNSKSNLLD TEHGDKVAEQ MNLASMDCLP DDLVQILYF LPTKEAISTS LLSKRWRTLY SLVHNLDD YIFWHHEDGY NQYFSNDIQK SFEEFMERTL ALLGGGHIKT FSLISDEIYR FDHVDSIRP WLYYNLQKDS LWQFGFPYKV FTSTKLVKLS LGTRLACPRI PQDTSPLVK VLLLEYIWFE DNQLSDVFLA ACPALEDLTI HHMFRPFLIS SKNLKLSVT INFSYYVDRS TILTLDTPNV VDLYYSDFPR PIAPHCHLDS LAKVELDLHS LEDDSRQVQN DADVKNLISE IRNVKTLHLT YSAVELMLSS KKRDKWVLP LERSPNLKT LVLSGLNRYT FGRHRFVGI QIPSSNKIKM LSIKQYQGS TELKHISHLL LKMECLELVK VYLATEMDDL KKMQLTEDVV KLPAASSKIV EKLKAFGAKA VLERVSFLD
Source	Mammalian cell
Target Names	At3g44090
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At3g44090
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.