



Recombinant Arabidopsis thaliana 50S ribosomal protein L1, chloroplastic (RPL1)

| | |
|--------------------------|---|
| Product Code | CSB-EP862929DOA |
| Abbreviation | RPL1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9LY66 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | IVS AVAAEADLDT EEDLEQTATA VLDPPKPKKG KAALVLKRDR TRSKRFLEIQ KLRETKKEYD VNTAISLLKQ TANTRFVESV EAHFRLNIDP KYNDQQLRAT VSLPKGTGQT VIVAVLAQGE KVDEAKSAGA DIVGSDDLIE QIKGGFMEFD KLIASPDMMV KVAGLGKILG PRGLMPNPKA GTVTANIPQA IEEFKKGKVE FRADKTGIVH IPFGKVNFTEDLLINFLAA VKSVETNPKPK GAKGVYWKSA HICSSMGPSI KLNIREMIDF KPPTAN |
| Source | E.coli |
| Target Names | RPL1 |
| Protein Names | Recommended name: 50S ribosomal protein L1, chloroplastic Alternative name(s): CL1 |
| Expression Region | 68-346 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |