



Recombinant Arabidopsis thaliana Adenosine kinase 2 (ADK2)

Product Code	CSB-MP862941DOA
Abbreviation	ADK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LZG0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASSSNYDGI LLGMGNPLLD ISAVVDDEFLL TKYDIKLNNA ILAEDKHLPM YDEMSSKFNV EYIAGGATQN SIKVAQWMLQ IPGATSYMGS IGKDKYGEAM KKDATAAGVN VHYYEDESAP TGTCGVCVVG GERSLIANLS AANCYKVDHL KKPENWALVE KAKFYIAGF FLTVSPESIQ LVSEHAAANN KVFTMNL SAP FICEFFKDVQ EKFLPYMDFV FGNETEARTF SRVHWETED VEQIAIKISQ LPKATGTYKR TTVITQGADP VVVAEDGKVK KYPVIPLPKE KLVD TNGAGD AFVGGFMSQL VKEKSIEECV KAGCYASN VV IQRSGCTYPE KP DFN
Source	Mammalian cell
Target Names	ADK2
Protein Names	Recommended name: Adenosine kinase 2 Short name= AK 2 EC= 2.7.1.20 Alternative name(s): Adenosine 5'-phosphotransferase 2
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.