



Recombinant Arabidopsis thaliana Putative two-component response regulator ARR19 (ARR19)

Product Code	CSB-EP863002DOA-B
Abbreviation	ARR19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M9B9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLVGKISGYE DNTRSLERET SEITSLLSQF PGNTNVLVVD TNFTLLNMK QIMKQYAYQV SIETDAEKAL AFLTSCKHEI NIVIWDFHMP GIDGLQALKS ITSKLDLPVV IMSDDNQTES VMKATFYGAC DYVVKPVKEE VMANIWQHIV RKRLIFKPDV APPKPRMTWS EVFQPVQSHL VPTDGLDRDH FDSITINGGN GIQNMEKKQG KKPRKPRMTW TEELHQKFLE AIEIIGGIEK ANPKVLVECL QEMRIEGITR SNVASHLQKH RINLEENQIP QQTQGNQWAT AYGTLAPSLQ GSDNVNTTIP SYLMNGPATL NQIQQNQYQN GFLTMNNNQI ITNPPPPLPY LDHHHQQQHQQ SSPQFNLYLMN NEELLQASGL SATDLELTYP SLPYDPQEYL INGYNYN
Source	E.coli
Target Names	ARR19
Protein Names	Recommended name: Putative two-component response regulator ARR19
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.