



Recombinant *Oryza sativa* subsp. *japonica* Coatomer subunit epsilon-1 (COPE1)

Product Code	CSB-MP863015OFG
Abbreviation	COPE1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9MAX6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASPDLLFNL RNLFYLGAYQ AAINNSDVPG LDADAAAERD AIVFRSYVAL GSYQLVISEI DSSAATSLQA VKLLALYLSG DKESAIVSLK EWLSDSAVGS NPVLRLIAGI IFMHEQDYTE ALKHTHSGGT LDLHALNVQI FIKMHRSDYA EKQLKIMQQI DEDHTLTQLA NAWLDIAVGG SKIREAYLIF QDFAEKYPMT GMVLNGKAVC CMHMGSFDEA ETLLEALNK DAKDPETLAN LIVCNLHLGK PSSRYLSQLK LSHPDHVLVK RAVSAEDNFE RALQAVA
Source	Mammalian cell
Target Names	COPE1
Protein Names	Recommended name: Coatomer subunit epsilon-1 Alternative name(s): Epsilon-coat protein 1 Short name= Epsilon-COP 1 Epsilon1-COP
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.