



Recombinant Cat Interferon beta (IFNB1)

Product Code	CSB-YP863088CA
Abbreviation	IFNB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N2J0
Product Type	Recombinant Protein
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Purity	≥85% (SDS-PAGE)
Sequence	VSYKLLGFQ LRSSSLECQE LLVNLNRTSK YCLKDRMNFE VP EEIKKSQR FQKEEAILVI NEMFQKIFNI FSRSTSSTGW NETTVENLLA TLHWQKEHLE TILEEIMEEE NFTWDNTLL NLKKYYLRIV RYLKAKEYSV CAWTVVHAEI LRNFFFLERL TDYLQN
Source	Yeast
Target Names	IFNB1
Protein Names	Recommended name: Interferon beta Short name= IFN-beta
Expression Region	22-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.