



Recombinant 60S ribosomal protein L10a (rpl-10a)

Product Code	CSB-YP863094CXY
Abbreviation	rpl-10a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N4I4
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSKVSRESLN EAIAEVLKGS SEKPRKFRET IELQIGLKNY DPQKDKRFSG SIRLKHIPRP NMKVCVFGDQ HHLDEAAAGD IPSMSADDLK KLNKQKKLIK KLAksYDAFI ASESliKQIP RILGPGLNKA GKFPSVvTHG ESLQSKSDEI RATVkfQMkk VLCLsvAVGH VGLTQEELVS NISLSINFLV SLLKKNWQNV RSLNIKSTMG KPQRVY
Source	Yeast
Target Names	rpl-1
Protein Names	Recommended name: 60S ribosomal protein L10a
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.