



Recombinant Human Dysbindin domain-containing protein 2 (DBNDD2)

Product Code	CSB-MP863102HU
Abbreviation	DBNDD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQY9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGAGNFLTAL EVPVAALAGA ASDRRASCER VSPPPPLPHF RLPPLPRSRL PGPVSRRPEPG APLLGCWLQW GAPSPGPLCL LFRLCSTCF APLPAGADMD PNPRAALERQ QLRLRERQKF FEDILQPETE FVFPLSHLHL ESQRPIIGSI SSMEVNVDTL EQVELIDLGD PDAADVFLPC EDPPPTPQSS GMDNHLEELS LPVPTSDRTT SRTSSSSSSD SSTNLHSPNP SDDGADTPLA QSDEEEERGD GGAEPGACS
Source	Mammalian cell
Target Names	DBNDD2
Protein Names	Recommended name: Dysbindin domain-containing protein 2 Short name= Casein kinase-1 binding protein Alternative name(s): CK1BP HSMNP1
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.