



# Recombinant Human BUD13 homolog (BUD13)

<b>Product Code</b>	CSB-MP863104HU
<b>Abbreviation</b>	BUD13
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BRD0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MAAAPLSKA EYLKRYLSGA DAGVDRGSES GRKRRKKRPK            PGGAGGKGMR IVDDVSWTA ISTTKLEKEE EEDDGDLPVV AEFVDERPEE            VKQMEAFRSS AKWKLLGGHN EDLPSNRHFR HDTDPSSPRR            VRHGTPDPSP RKDRHDTDP SPRRARHDTP DPSPLRGARH DSDTSPPRRI            RHDSSDTSP RRRARHDSPDP SPPRRPQHNS SGASPRRVRH            DSPDPSPRR ARHGSSDISS PRRVHNNSPD TSRRTLGSDD TQQLRRARHD            SPDLAPNVTY SLPRTKSGKA PERASSKTSP HWKESGASHL SFPKNSKYEY            DPDISPPRKK QAKSHFGDKK QLDSKGCQK ATDSDLSSPR HKQSPGHQDS            DSDLSPPRNR PRHRSSDSL SPPRRRQRTK SSDSDLSPPR RSQPPGKAA            HMYSGAKTGL VLTDIQREQQ ELKEQDQETM AFEAEFQYAE TVFRDKSGRK            RNLKLERLEQ RRKAEKDSER DELYAQWGKG LAQSRQQQQN            VEDAMKEMQK PLARYIDDED LDRMLREQER EGDPMANFIK KNKAKENKNK            KVRPRYSGPA PPPNRFNIWP GYRWDGVDRS NGFEQKRFAR LASKKAVEEL            AYKWSVEDM         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BUD13
<b>Protein Names</b>	Recommended name: BUD13 homolog
<b>Expression Region</b>	1-619
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.