



Recombinant Human BUD13 homolog (BUD13)

Product Code	CSB-YP863104HU
Abbreviation	BUD13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BRD0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAAAPPLSKA EYLKRYLSGA DAGVDRGSES GRKRRKKRPK PGGAGGKGMR IVDDVSWTA ISTTKLEKEE EEDDGDLPVV AEFVDERPEE VKQMEAFRSS AKWKLLGGHN EDLPSNRHFR HDTDPSSPRR VRHGTPDPSP RKDRHDTDP SPRRARHDTP DPSPLRGARH DSDTSPPRRI RHDSSDTSP RRRARHDSPDP SPPRRPQHNS SGASPRRVRH DSPDPSPRR ARHGSSDISS PRRVHNNSPD TSRRTLGSDD TQQLRRARHD SPDLAPNVTY SLPRTKSGKA PERASSKTSP HWKESGASHL SFPKNSKYEY DPDISPPRKK QAKSHFGDKK QLDSKGDQK ATDSDLSSPR HKQSPGHQDS DSDLSPPRNR PRHRSSDSL SPPRRRQRTK SSDSDLSPPR RSQPPGKAA HMYSGAKTGL VLTDIQREQQ ELKEQDQETM AFEAEFQYAE TVFRDKSGRK RNLKLERLEQ RRKAEKDSER DELYAQWGKG LAQSRQQQQN VEDAMKEMQK PLARYIDDED LDRMLREQER EGDPMANFIK KNKAKENKNK KVRPRYSGPA PPPNRFNIWP GYRWDGVDRS NGFEQKRFAR LASKKAVEEL AYKWSVEDM</p>
Source	Yeast
Target Names	BUD13
Protein Names	Recommended name: BUD13 homolog
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.