



Recombinant Human WD repeat domain-containing protein 83 (WDR83)

Product Code	CSB-BP863110HU
Abbreviation	WDR83
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BRX9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAFPEPKPRP PELPQKRLKT LDCGQGAVRA VRFNVDGNYC LTCGSDKTLK LWNPLRGTLT RTYSGHGYEV LDAAGSFDNS SLCSGGGDKA VVLWDVASGQ VVRKFRGHAG KVNTVQFNEE ATVILSGSID SSIRCWDCRS RRPEPVQTLT EARDGVSSVK VSDHEILAGS VDGRVRRYDL RMGQLFSDYV GSPITCTCFS RDGQCTLVSS LDSTLRLLDK DTGELLGEYK GHKNQEYKLD CCLSERDTHV VSCSEDGKVF FWDLVEGALA LALPVGSGVV QSLAYHPTEP CLLTAMGGSV QCWREEAYEA EDGAG
Source	Baculovirus
Target Names	WDR83
Protein Names	Recommended name: WD repeat domain-containing protein 83 Alternative name(s): Mitogen-activated protein kinase organizer 1 Short name= MAPK organizer 1
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.