



Recombinant Human Polycomb group RING finger protein 1 (PCGF1)

Product Code	CSB-YP863114HU
Abbreviation	PCGF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSM1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASPQGGQIA IAMRLRNQLQ SVYKMDPLRN EEEVRVKIKD LNEHIVCCLC AGYFVDATTI TECLHTFCKS CIVKYLQTSK YCPMCNIKIH ETQPLLNLKL DRVMQDIVYK LVPGLQDSEE KRIREFYQSR GLDRVTQPTG EEPALSNLGL PFSSFDHKA HYYRYDEQLN LCLERLSSGK DKNKSVLQNK YVRCSVRAEV RHLRRVLCHR LMLNPQHVQL LFDNEVLPDH MTMKQIWLRS WFGKPSPLLL QYSVKEKRR
Source	Yeast
Target Names	PCGF1
Protein Names	Recommended name: Polycomb group RING finger protein 1 Alternative name(s): Nervous system Polycomb-1 Short name= NSPc1 RING finger protein 68
Expression Region	2-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.