



Recombinant Human HAUS augmin-like complex subunit 8 (HAUS8)

Product Code	CSB-YP863116HU
Abbreviation	HAUS8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT25
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ADSSGRGAG KPATGPTNSS SAKKKDKRVQ GGRVIESRYL QYEKTTQKA PAGDGSQTRG KMSEGGRKSS LLQKSKADSS GVGKGDQLQST LLEGHTAPP DLDLSAINDK SIVKKTPQLA KTISKKPEST SFSAPRKKSP DLSEAMEMME SQTLLLTLLS VKMENNLAEF ERRAEKNNLI MCKEKEKLQK KAHELKRRL LSQRKRELAD VLDAQIEMLS PFEAVATRFK EQYRTFATAL DTTRHELPVR SIHLEGDGQQ LLDALQHELV TTQRLLGELD VGDSEENVQV LDLLSELKDV TAKKDLELRR SFAQVLELSA EASKEAALAN QEVWEETQGM APPSRWYFNQ DSACRESGGA PKNTPLSEDD NPGASSAPAQ ATFISPSEDF SSSSQAEVPP SLSRSGRDLS
Source	Yeast
Target Names	HAUS8
Protein Names	Recommended name: HAUS augmin-like complex subunit 8 Alternative name(s): HEC1/NDC80-interacting centrosome-associated protein 1 Sarcoma antigen NY-SAR-48
Expression Region	2-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.