



Recombinant Human tRNA-specific adenosine deaminase 1 (ADAT1)

Product Code	CSB-EP863123HU
Abbreviation	ADAT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BUB4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MWTADEIAQL CYEHYGIRLP KKGKPEPNHE WTLAAVVKI QSPADKACDT PDKPVQVTKE VVSMGTGTKC IGQSKMRKNG DILNDSHAEV IARRSFQRYL LHQLQLAATL KEDSIFVPGT QKGVWKLRRD LIFVFFSST PCGDASIIPM LEFEDQPCCP VFRNWAHNS VEASSNLEAP GNERKCEDPD SPVTKKMRLE PGTAAREVTN GAAHHQSFGK QKSGPISPGI HSCDLTVEGL ATVTRIAPGS AKVIDVYRTG AKCVPGEAGD SGKPGA AFHQ VGLLRVKPGR GDRTRSMSCS DKMARWNVLG CQGALLMHLL EEPIYLSAVV IGKCPYSQEA MQRALIGRCQ NVSALPKGFG VQELKILQSD LLFEQSRSAV QAKRADSPGR LVPCGAAISW SAVPEQPLDV TANGFPQGTT KKTIGSLQAR SQISKVELFR SFQKLLSRIA RDKWPHSLRV QKLDTYQEYK EAASSYQEA WSTLRKQVFGS WIRNPPDYHQ FK
Source	E.coli
Target Names	ADAT1
Protein Names	Recommended name: tRNA-specific adenosine deaminase 1 Short name= hADAT1 EC= 3.5.4.- Alternative name(s): tRNA-specific adenosine-37 deaminase
Expression Region	1-502
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the ADAR (adenosine deaminase acting on RNA) family. Using site-specific adenosine modification, proteins encoded by these genes participate in the pre-mRNA editing of nuclear transcripts. This protein, tRNA-specific adenosine deaminase 1, is responsible for the deamination of adenosine 37 to inosine in eukaryotic tRNA.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.