



Recombinant Human Probable tubulin polyglutamylase TTLL2 (TTLL2)

Product Code	CSB-MP863137HU
Abbreviation	TTLL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWW7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRGRDLCSSST QSQALGSLRT TTPAFTLNIP SEANHTEQPP AGLGARLQEA GVSIPRRGR PTPTLEK KKK PHLMAEDEPS GALLKPLVFR VDETTPAVVQ SVLLERGWNK FDKQE QNAED WNL YWRTSSF RMTEHNSVKP WQQLNHHPGT TKLTRKDCLA KHLKHMRRMY GTSLYQFIPL TFVMPNDYTK FVAEYFQERQ MLGTKHSYWI CKPAELSRGR GILIFSDFKD FIFDDMYIVQ KYISNPLLIG RYKCDLRIYV CVTGFKPLTI VVYQEGLVRF ATEKFDLSNL QNNYAHLTNS SINKSGASYE KIKEVIGHGC KWTL SRFFSY LRSWDVDDLL LWKKIHRMVI LTI LAIPSV PFAANCFELF GFDILIDDNL KPWLLEVNYS PALTLDCSTD VLVKRKLVHD IIDLIYLNGL RNEGREASNA THGNSNIDAA KSDRGGLDAP DCLPYDSL SF TSRMYNEDDS VVEKAVSVRP EAAPASQLEG EMSGQDFHLS TREMPQSKPK LRSRHTPHKT LMPYASLFQS HSCKTKTSPC VLSDRGKAPD PQAGNFVLVF PFNEATLGAS RNGLNVKRII QELQKLMNKQ HS
Source	Mammalian cell
Target Names	TTLL2
Protein Names	Recommended name: Probable tubulin polyglutamylase TTLL2 EC= 6.-.-.- Alternative name(s): Testis-specific protein NYD-TSPG Tubulin--tyrosine ligase-like protein 2
Expression Region	1-592
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.