



Recombinant Human Rab11 family-interacting protein 5 (RAB11FIP5)

Product Code	CSB-MP863141HU
Abbreviation	RAB11FIP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXF6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALVRGAEPA AGPSRWLPTH VQVTVLRARG LRGKSSGAGS TSDAYTVIQV GREKYSTSVV EKTHGCPEWR EECSEFELPPG ALDGLLRAQE ADAGPAPWAA SSAAACELVL TTMHRSLIGV DKFLGQATVA LDEVFGAGRA QHTQWYKLHS KPGKKEKERG EIEVTIQFTR NNLSASMF DL SMKDKPRSPF SKIRDKMKGK KKYDLESASA ILPSSAIEDP DLGSLGKMGK AKGFFLRNKL RKSSLTQSNT SLGSDSTLSS ASGSLAYQGP GAELLTRSPS RSSWLSTEGG RDSAQSPKLF THKRTYSDEA NQMRVAPPRA LLDLQGH LDA ASRSSLCVNG SHIYNEEPQG PVRHRSSISG SLPSSGSLQA VSSRFSEEGP RSTDDTWPRG SRSNSSSEAV LGQEELSAQA KVLAPGASHP GEEEGARLPE GKPVQVATPI VASSEAVA EK EGARKEERKP RMGLFH HH HQ GLSRSELGRR SSLGEKGGPI LGASPHHSSS GEEKAKSSWF GLREAKDPTQ KPSPHPVKPL SAAPVEGSPD RKQSRSSLSI ALSSGLEK LK TVTSGSIQPV TQAPQAGQMV DTKRLKDSAV LDQSAKY YHL THDELISLLL QRERELSQRD EHVQELESYI DRLLVRIMET SPTLLQIPPG PPK
Source	Mammalian cell
Target Names	RAB11FIP5
Protein Names	Recommended name: Rab11 family-interacting protein 5 Short name= Rab11-FIP5 Alternative name(s): Gamma-SNAP-associated factor 1 Short name= Gaf-1 Phosphoprotein pp75 Rab11-interacting protein Rip11
Expression Region	1-653
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.