



Recombinant Human 28S ribosomal protein S26, mitochondrial (MRPS26)

Product Code	CSB-MP863157HU
Abbreviation	MRPS26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYN8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	KTR HDPLAKSKIE RVNMPPAVDP AEFVLMERY QHYRQTVRAL RMEFVSEVQR KVHEARAGVL AERKALKDAA EHRELMAWNQ AENRRRLHELRL IARLRQEERE QEQRQALEQA RKAEEVQAWA QRKEREVLQL QEEVKNFITR ENLEARVEAA LDSRKNYNWA ITREGLVVRP QRRDS
Source	Mammalian cell
Target Names	MRPS26
Protein Names	Recommended name: 28S ribosomal protein S26, mitochondrial Short name= MRP-S26 Short name= S26mt Alternative name(s): 28S ribosomal protein S13, mitochondrial Short name= MRP-S13 Short name= S13mt
Expression Region	28-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. This gene lies adjacent to and downstream of the gonadotropin-releasing hormone precursor gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.