



Recombinant Human tRNA pseudouridine (38/39) synthase (PUS3)

Product Code	CSB-EP863167HU
Abbreviation	PUS3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BZE2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AYNDTDRNQ TEKLLKRVRE LEQEVQRLKK EQAKNKEDSN IRENSAGAGK TKRAFDFSAH GRRHVALRIA YMGWGYQGFA SQENTNNTIE EKLFEALTKT RLVESRQTSN YHRCGRTDKG VSAFGQVISL DLRSQFPRGR DSEDFNVKEE ANAAEEIERY THILNRVLPP DIRILAWAPV EPSFSARFSC LERTYRYFFP RADLDIVTMD YAAQKYVGTH DFRNLCKMDV ANGVINFQRT ILSAQVQLVG QSPGEGRWQE PFQLCQFEVT GQAFLYHQVR CMMAILFLIG QGMEKPEIID ELLNIEKNPQ KPQYSMAVEF PLVLYDCKFE NVKWIYDQEA QEFNITHLQQ LWANHAVKTH MLYSMLQGLD TVPVP CGIGP KMDGMTEWGN VKPSVIKQTS AFVEGVKMRT YKPLMDRPKC QGLESRIQHF VRRGRIEHPH LFHEEETKAK RDCNDTLEEE NTNLETPTKR VCV DTEIKSI I
Source	E.coli
Target Names	PUS3
Protein Names	Recommended name: tRNA pseudouridine(38/39) synthase EC= 5.4.99.45 Alternative name(s): tRNA pseudouridine synthase 3 tRNA pseudouridylate synthase 3 tRNA-uridine isomerase 3
Expression Region	2-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.