



# Recombinant Human tRNA pseudouridine (38/39) synthase (PUS3)

<b>Product Code</b>	CSB-MP863167HU
<b>Abbreviation</b>	PUS3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BZE2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AYNDTDRNQ TEKLLKRVRE LEQEVQRLKK EQAKNKEDSN IRENSAGAGK TKRAFDFSAH GRRHVALRIA YMGWGYQGFA SQENTNNTIE EKLFEALTKT RLVESRQTSN YHRCGRTDKG VSAFGQVISL DLRSQFPRGR DSEDFNVKEE ANAAEEIERY THILNRVLPP DIRILAWAPV EPSFSARFSC LERTYRYFFP RADLDIVTMD YAAQKYVGTH DFRNLCKMDV ANGVINFQRT ILSAQVQLVG QSPGEGRWQE PFQLCQFEVT GQAFLYHQVR CMMAILFLIG QGMEKPEIID ELLNIEKNPQ KPQYSMAVEF PLVLYDCKFE NVKWIYDQEA QEFNITHLQQ LWANHAVKTH MLYSMLQGLD TVPVP CGIGP KMDGMTEWGN VKPSVIKQTS AFVEGVKMRT YKPLMDRPKC QGLESRIQHF VRRGRIEHPH LFHEEETKAK RDCNDTLEEE NTNLETPTKR VCV DTEIKSI I
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PUS3
<b>Protein Names</b>	Recommended name: tRNA pseudouridine(38/39) synthase EC= 5.4.99.45 Alternative name(s): tRNA pseudouridine synthase 3 tRNA pseudouridylate synthase 3 tRNA-uridine isomerase 3
<b>Expression Region</b>	2-481
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.