



Recombinant Human Testis-specific gene 10 protein (TSGA10)

Product Code	CSB-BP863172HU
Abbreviation	TSGA10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BZW7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMRSRSKSPR RPSPTARGAN CDVELLKTTT RDREELKCML EKYERHLAEI QGNVKVLKSE RDKIFLLYEQ AQEEITRLRR EMMKSCKSPK STTAHAILRR VETERDVAFT DLRRMTTERD SLRERLKIAQ ETAFNEKAHL EQRIEELECT VHNLDDERME QMSNMTLMKE TISTVEKEMK SLARKAMDTE SELGRQKAEN NSLRLLYENT EKDLSDTQRH LAKKKYELQL TQEKIMCLDE KIDNFTRQNI AQREEISILG GTLNDLAKEK ECLQACLDKK SENIASLGES LAMKEKTISG MKNIIAEMEQ ASRQCTEALI VCEQDVSRMR RQLDETNDL AQIARERDIL AHDNDNLQEQ FAKAKQENQA LSKKLNDTHN ELNDIKQKVQ DTNLEVNKLK NILKSEESN RQMMEQLRKA NEDAENWENK ARQSEADNNT LKLELITAE EGNRLKEKVD SLNREVEQHL NAERSYKSQI STLHKS SVVKM EEELQKVQFE KVSALADLSS TRELCIKLDS SKELLNRQLV AKDQEIE MRE NELDSA HSEI ELLRSQMANE RISMQNLEAL LVANRDKEYQ SQIALQEKES EIQLLKEHLC LAENKMAIQS RDVAQFRNVV TQLEADLDIT KRQLGTERFE RERAVQELRR QNYSSNAYHM SSTMKPNTKC HSPERAHHRS PDRGLDRSLE ENLCYRDF
Source	Baculovirus
Target Names	TSGA10
Protein Names	Recommended name: Testis-specific gene 10 protein Alternative name(s): Testis development protein NYD-SP7
Expression Region	1-698
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.