



Recombinant Schizosaccharomyces pombe UPF0613 protein PB24D3.06c (SPAPB24D3.06c)

Product Code	CSB-MP863186SXV
Abbreviation	SPAPB24D3.06c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0Y8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSASHPGILH EYTERLVAFE LVGKSNVTLT SNVTRSLLLF VGGLGDGLLT VPYVQELVNP LDEIGWSIVQ VQTQSSYIGW GTGSLKRDDE DLHKAVDYFL HIGGADFSTR KIVLMGHSTG SQNVLYYLTQ SILPNYLIAG IAQAPVSDRE AAYQFNGKEK TKELVDWVKA EYLDKGLGND VLPRSKVENF FGEVPTSANR CIDLTDVRGN DDFSSDLA DDFAKTFGNL KEISGSTAHS QLILLMSERD EFVSPSTDKA QLLNRFRESI RPTTSSSSLS GIIPGATHNV GPKSSPEALK WLINQLITAL RSFIDQ
Source	Mammalian cell
Target Names	SPAPB24D3.06c
Protein Names	Recommended name: UPF0613 protein PB24D3.06c
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.