



Recombinant Arabidopsis thaliana UPF0725 protein At1g27860 (At1g27860)

Product Code	CSB-YP863218DOA
Abbreviation	At1g27860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6N2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPTEDRVCFT WEEYARHVRE YWRGVAESDG FDSEDIRSPV VLTGLFHYDC QSGRHYPEPL LGTSFELAAL LKFNKTMNLT SSYYITLLAH DPSLEKTFQV RVDEREYGSL DLTVAIARPK KDQNEVSDGV FQGPLPDWPS EDALRHDRNR FYELEISEWQ ATDWISLYLK LLILATDRGM FVQTGLPQVQ ILEVVIETEE ENEKPPDKRL NARRAHVYIT FTGLPKSPRL VEIGEHVERK AIIRRVIDAS GYLTLGKFW SGKDTDQQRS MTLPSGEQAE SSKKRPRLS
Source	Yeast
Target Names	At1g27860
Protein Names	Recommended name: UPF0725 protein At1g27860
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.