



Recombinant Arabidopsis thaliana Transcription factor TCP22 (TCP22)

Product Code	CSB-BP863223DOA
Abbreviation	TCP22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C7G4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNQNSSVAEA TLQLNSGEKP SPGSIPFISS GQHGNISTSA TSSTSTSSGS ALAVVKS AVK KPTKDRHTKV DGRGRRIRMP AMCAARVFQL TRELGHKSDG ETIEWLLQQA EPAIIASTGT GTIPANFSTL NASLRSGGGS TLFSQASKSS SSPLSFHSTG MSLYEDNNGT NGSSVDPSRK LLNSAANA AV FGFHHQMYPP IMSTERNPNT LVKPYREDYF KEPSSAAEPS ESSQKASQFQ EQELAQQGRGT ANVVPQPMWA VAPGTTNGGS AFWMLPMSGSG GREQMQQQP GHQMWAFNPG NYPVGTGRVV TAPMGSMMLG GQQLGLGVAE GNMAAAMRGS RGDGLAMTLD QHQHQLQHQE PNQSQASENG GDDKK
Source	Baculovirus
Target Names	TCP22
Protein Names	Recommended name: Transcription factor TCP22
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.