



Recombinant Arabidopsis thaliana Putative zinc finger CCCH domain-containing protein 10 (At1g29600/At1g29610)

Product Code	CSB-EP863224DOA-B
Abbreviation	At1g29600/At1g29610
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C7P1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MANVSFTFDS QEQNKELRPS RLVGRSDERN GRRQGTEQKK KEKEDSGFRV SEEGEIPLQQ GSNTDSYARD DLTQQQRVSD VVMQRRSHET ESRLWQRART QDRRGSESRM MFDGRTQWSH APVLSSAYPV RPGEDNCLFY MKNHLCWGS ECCYNHPPLQ EIPCRIGKKL DCKAGACKRG SNCPFNHPKE RDGDSLMPQ GRTPDLRRND SGRRYNTESR SWPENKEKEV GQFRDHQDSK EDAQEVLLQQ RPRDVEMRKR SRSPDFRAKT ETKEHREAEE RSSRESATAT GKVSGKENEL RPTYVYLLMM PENHNVDVQQ TLQEQRKRNI KKARIEARLK LDQIRPTVVL ENSRATKMMI EWGFTDFRKL QGSNLGNKKW ICWVRKHCC
Source	E.coli
Target Names	At1g29600/At1g29610
Protein Names	Recommended name: Putative zinc finger CCCH domain-containing protein 10 Short name= AtC3H10
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.