



# Recombinant Arabidopsis thaliana Putative zinc finger CCCH domain-containing protein 10 (At1g29600/At1g29610)

<b>Product Code</b>	CSB-MP863224DOA
<b>Abbreviation</b>	At1g29600/At1g29610
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C7P1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MANVSFTFDS QEQNKELRPS RLVGRSDERN GRRQGTEQKK KEKEDSGFRV SEEGEIPLQQ GSNTDSYARD DLTQQQRVSD VVMQRRSHET ESRLWQRART QDRRGSESRM MFDGRTQWSH APVLSSAYPV RPGEDNCLFY MKNHLCWGS ECCYNHPPLQ EIPCRIGKKL DCKAGACKRG SNCPFNHPKE RDGDSLMPQ GRTPDLRRND SGRRYNTESR SWPENKEKEV GQFRDHQDSK EDAQEVLLQQ RPRDVEMRKR SRSPDFRAKT ETKEHREAEE RSSRESATAT GKVSGKENEL RPTYVYLLMM PENHNVDVQQ TLQEQRKRNI KKARIEARLK LDQIRPTVVL ENSRATKMMI EWGFTDFRKL QGSNLGNKKW ICWVRKHCC
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At1g29600/At1g29610
<b>Protein Names</b>	Recommended name: Putative zinc finger CCCH domain-containing protein 10 Short name= AtC3H10
<b>Expression Region</b>	1-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.