



# Recombinant Arabidopsis thaliana Probable sucrose-phosphatase 1 (SPP1)

<b>Product Code</b>	CSB-MP863231DOA
<b>Abbreviation</b>	SPP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C8J4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MERLTSPRL MIVSDLDETM VDHKDPENL ALLRFNSLWE DAYRHDSLLV FSTGRAQTMV KCLRKEKPLL TPDVIITSVG TEIAYGNSMV PDDNWVEILN KKWDRGIVQE ETSKFPELTL QGETEQRPHK LSFNIDKSKV KAVTKELSPR LEKRGVDVKF IFSGGNAFDV LAKGGGKGQA LAYLLKLLKT EGKLPINTLA CGDSGNDTEL FTIPNVYGVV VSNAQEELLE WYAENAKDNA NIIHASERCA GGITQAIGHF KLGPNLSPRD VPDFLECKAD NVNPGHEVVK FFLFYERWRR GEVENCTTYT SSLKASCHPS GVFVHPSGAE KSLRDTIDEL GKYYGDKKGGK KFRVWTDQVL ATDTPGTWI VKLDKWEQTG DERKCCTTTV KFTSKEGEGF VWEHVQQIWS EETEIKDDSN WII
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPP1
<b>Protein Names</b>	Recommended name: Probable sucrose-phosphatase 1 Short name= AtSPP1 EC= 3.1.3.24
<b>Expression Region</b>	1-423
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.